



THE EVERFLO™ OXYGEN CONCENTRATOR
HAS BEEN DESIGNED TO MAKE A BIG IMPRESSION
WITH ITS SMALLER SIZE AND LIGHTER WEIGHT

THE SPACE-SAVING EVERFLO OXYGEN CONCENTRATOR



EverFlo™ Oxygen Concentrator

The smaller, lighter EverFlo Oxygen Concentrator from Respironics is making a big impression on providers and patients. Its sleek, ergonomic design makes EverFlo easier to store, move, and carry than conventional concentrators. And it is one of the quietest and most energy efficient devices available today. With no filter for patients to maintain, a durable metal oxygen outlet, and fewer moving parts than other concentrators, EverFlo is designed to deliver consistent, trouble-free operation.

Product Specifications

	Standard Version	Low Power Version	IKK Version
POWER	230V/50Hz	230V/50Hz	230V/50Hz
LITER FLOW	0.5-5 liters per min	0.5-5 liters per min	0.5-5 liters per min
SOUND LEVEL	43 dBA	43 dBA	< 40 dBA
OXYGEN PURITY (AT 5 LPM)	93% (+/- 3%)	90% (+/- 3%)	94% (+/- 2%) at 2 lpm 88% (+/- 4%) at 5 lpm
WEIGHT	14 kg (31 lbs)	14 kg (31 lbs)	14 kg (31 lbs)
DIMENSIONS (HxWxD)	580mm (22.8"), 380mm (15"), 240mm (9.5")	580mm (22.8"), 380mm (15"), 240mm (9.5")	580mm (22.8"), 380mm (15"), 240mm (9.5")
OPI (OXYGEN PERCENTAGE INDICATOR) ALARM LEVELS	Low Oxygen: 82%, Very Low Oxygen: 70%	Low Oxygen: 82%, Very Low Oxygen: 70%	Low Oxygen: 82%, Very Low Oxygen: 70%
REQUIRED MAINTENANCE	Compressor Intake Filter every 1-2 yrs. Compressor cup seals at 25,000 hrs.	Compressor Intake Filter every 1-2 yrs. Compressor cup seals at 25,000 hrs.	Compressor Intake Filter every 1-2 yrs. Compressor cup seals at 25,000 hrs.
POWER CONSUMPTION	295 Watts	270 Watts up to 3 lpm, 280 Watts at 5 lpm	295 Watts
PART NUMBERS	1020006 Europe 1020009 South Africa 1020011 Italy	1020008 UK, Ireland	1020007 Germany



Key features include:

- Lighter weight
- Lower power consumption
- Compact profile
- Lower sound level
- Five liter per minute oxygen delivery
- Two-year, low-maintenance filter



Respironics, Inc.: +1-724-387-4000
 Respironics Europe, Africa, Middle East: +33-1-47-52-30-00
 Respironics Deutschland: +49-81-52-930-60
 Respironics France: +33-2-51-89-36-00
 Respironics Italy: +39-03-62-63-43-1
 Respironics UK: +44-870-770-3434
 Respironics Asia Pacific: +852-3194-2280
www.respironics.com

Respironics is a registered trademark and EverFlo is a trademark of Respironics, Inc and its affiliates.
 ©2006 Respironics, Inc. and its affiliates. All rights reserved.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician.

BROUDY CY 11/22/06 MCI 4100822 PN 1039011