

prisma VENT40

Invasive and non-invasive ventilation.

With prisma VENT40 you'll have **mobility and comfort at every level:**

- Pressure range up to 40 hPa
- Use of leakage circuit system (15 mm / 22 mm)
- With autoST mode (autoEPAP, autoF) and target volume (100 – 2,000 ml)
- Includes helpful functions for COPD treatment: AirTrap Control, maximum setting for expiratory ramp and Trigger Lockout
- Integrated port for oxygen feed
- Three ventilation programs stored in device memory
- Complete alarm management
- Rechargeable battery with 12-hour capacity
- Easy to use thanks to intuitive menus, simple operation and quick access ...

... and it's **quiet as a whisper for the well-being and comfort of patients and partners.**

Our Software Solutions



- **prisma CLOUD**
The telemedical connection for experts and patients. Information updated daily on therapy, fully automated, from anywhere.
- **prismaTS**
Therapy software for experts. For detailed data analysis and remote settings.

