

Now also available in **XL**

Color-coded headgear – colored (teal or syrah) part on top

Soft, supple mask cushion in four sizes

Quick-release cord for rapid removal

Extremely durable material

Choice of manual threading or permanent assembly

Quiet and diffuse exhalation system

Ball-and-socket joint with 360° and 3D rotation for freedom of movement



Adjust the lower portion of the headgear the same way so that it is as loose as possible and as tight as necessary.



Position the headgear on the head so that the lower portion does not hit the ear lobes.



Adjust the forehead support as loosely as possible; it does not have to lie on the forehead.



You'll find our user video on YouTube.