

# JOYCEclinic Full Face

## No compromises.

With our **single-patient-use full-face mask** for acute non-invasive ventilation, you will experience the familiar quality from Löwenstein Medical.

JOYCEclinic Full Face allows **fast and simple mask fittings** and a **secure, comfortable fit** even at pressures of up to 50 hPa.

The headgear, which remains firmly in place even at high pressures, can be disconnected quickly and simply by pulling on the highly visible **quick-release cord**.

**Single-patient-use mask** prevents cross-contamination. In the hospital where cost pressure is high, you needn't make any compromises with patient safety.

- Prevention of cross-contamination
- Pressurestable mask cushion and headgear
- Quick-release cord

With four different, interchangeable elbow adapters, you'll find the right mask variant for every situation and device combination.

**JOYCEclinic Full Face** comes preassembled in three sizes in the variants **non-vented**, **non-vented with AAV** and **vented with AAV**. You can use the fitting template to make a good fit quickly and reliably. You don't even touch the mask, thanks to the two-chamber bag.

Available as a practical accessory, the **endoscope adapter** permits bronchoscopy without interrupting ventilation. The elbow adapters NV, NV+ AAV and Vented also are available as **replacement parts**.

All masks are available in sets of 10 without the two-chamber bag. Each set comes with Instructions for use and a fitting template.

- Three sizes
- Four interchangeable elbows
- Fitting template
- Two-chamber bag
- Single mask or set of 10 masks

