

prisma SOFT plus / prisma SMART plus prisma SOFT max / prisma SMART max

CPAP and APAP therapy: Flexible. Secure. Reliable.

With prisma SMART and prisma SOFT in the **plus and max variants**, we offer modern, reliable CPAP and APAP therapy featuring integrated connectivity, precise event detection and a high level of patient comfort.

The **plus** stands for integrated Bluetooth for data transmission to the prisma APP and thus for modern patient self-management. Equipped with an integrated modem, the **max** variants also offer the complete telemedical prisma CLOUD connection.

The most important features at a glance

- Two dynamic options for pressure adjustment in APAP mode
- Deep sleep indicator prisma RECOVER for evaluation of sleep quality
- Several convenience functions such as autoSTART, softSTART, pressure relief softPAP
- Detection of periodic breathing, RERA, snoring, hypopnea, apnea and flow limitations
- Reliable differentiation of obstructive apnea from central apnea, thanks to FOT (Forced Oscillatory Technique) ...

... and it's **quiet as a whisper for the well-being and comfort of patients and partners.**

Our Software Solutions



- **prisma CLOUD**
The telemedical connection for experts and patients. Information updated daily on therapy, fully automated, from anywhere.
- **prisma APP**
The prisma APP supports and motivates patients to adhere to treatment.
- **prismaTS**
Therapy software for experts. For detailed data analysis and remote settings.

