

BiPAP S/T



Noninvasive home ventilator

Physical

Size 18 cm x 14 cm x 10 cm (7" L x 5.5" W x 4" H)

Weight Approximately 1.36 kg (3 lbs.)

Therapy parameters

Modes CPAP, S, ST

Bi-Flex Available on S mode

Timed inspiration 0.5 to 3.0 seconds

IPAP 4 to 25 cm H₂O

EPAP 4 to 25 cm H₂O

CPAP 4 to 20 cm H₂O

Ramp 0 to 45 minutes

Breath rate 0 to 30 BPM

Patient Disconnect, Apnea, Low Minute Ventilation

Environmental	Operating	Storage
Temperature	5° C to 35° C (41° F to 95° F)	-20° C to 60° C (-4° F to 140° F)
Relative humidity	15 to 95% (non-condensing)	15 to 95% (non-condensing)
Atmospheric pressure	101 kPa to 77 kPa (0-2286 m / 0-7500 ft)	N/A

Sound level

< 30 dBA typical at 10 cm H₂O pressure

Electrical

AC voltage source 100 to 240 VAC, 50/60 Hz, 2.1 A

DC power 12 VDC, 5.0 A

Displayed parameters

Patient pressure, leak, tidal volume, minute ventilation, and respiratory rate

Filters

Foam and optional ultra-fine
